



Basic scATAC-seq workflow

Main Steps for scATACseq Workflow




Step1. Data processing



- Sample preparation involves tissue dissociation and nuclei isolation.
- **Library construction:** Tn5 transposition and barcoding.
- **Alignment:** Mapping reads to the reference genome.


Step1. Data processing



	Cell1	Cell2	...	CellN
Peak1	3	1	.	2
Peak2	2	0	.	1
Peak3	1	1	.	0
...
PeakM	5	1	.	0

- **Peak Calling:** Identifying accessible chromatin regions
- **Fragment Counting:** Quantifying accessibility per region

Step2. Quality control




- Why QC is crucial in snATAC-seq: Ensures high-quality and biologically relevant data.
- Common QC Metrics: Read depth, TSS enrichment, fragment size distribution, and nucleosome positioning.
- Doublet detection

Step3. Matrix formation & Normalization

	Cell1	Cell2	...	CellN		Cell1	Cell2	...	CellN
Peak1	3	1	.	2		0.1	0.2	.	0.3
Peak2	2	0	.	1		0.05	0	.	0.2
Peak3	1	1	.	0		0.02	0.2	.	0
...
PeakM	5	1	.	0		0.2	0.2	.	0


- Peak x Cell Matrix
- Cell-level Normalization: Adjusting for variations in sequencing depth across individual cells
- Peak-level Normalization: Scaling peak signals to assign higher values to rare peaks, enhancing the detection of less abundant regulatory regions

Step4. Cell clustering and annotation



- Cell Clustering: clustering after dimension reduction
- Cluster-level Annotation: using well-established cell type-specific markers
- Downstream Analysis: cluster analysis, differential accessibility analysis, and integration with multi-omics datasets

Step5. Downstream Analysis




Hands-On Example Dataset

Assay: Single Cell Multiome (ATAC + Gene Expression)

Sample Type: Peripheral Blood Mononuclear Cells (PBMCs)

Donors: Two healthy individuals

10x GENOMICS



- a healthy female donor aged 25

<https://www.10xgenomics.com/datasets/pbmc-from-a-healthy-donor-no-cell-sorting-10-k-1-standard-2-0-0>

10xGENOMICS

ProductsResourcesSupport HubCompany

Order StatusStoreSearch

[< All datasets](#)

PBMC from a Healthy Donor - No Cell Sorting (10k)


Epi Multiome ATAC + Gene Expression dataset analyzed using Cell Ranger ARC 2.0.0

Mapped
95.7%

Assess data quality

View summary metrics to assess data quality and more.

View summary



Visualize and explore data

Discover differentially expressed genes, visualize your favorite genes, and explore your data with our visualization software.

Explore data

Dataset overview

Output and supplemental files

Input files

[Learn about Chromium analysis](#)

Batch download

Input files	File type	Size	md5sum
Sequencing data (FASTQ)	TAR	114 GB	dc5b7f0c6bc3ace4b9c02780bdab32f5
Library (CSV)	CSV	180 B	90778b6daa76412f77d18cf9d51edf9a

- a healthy male donor aged 30-35

<https://www.10xgenomics.com/datasets/10-k-human-pbm-cs-multiome-v-1-0-chromium-controller-1-standard-2-0-0>

10X GENOMICS

ProductsResourcesSupport HubCompany

Order StatusStoreSearch

[< All datasets](#)

10k Human PBMCs, Multiome v1.0, Chromium Controller


Epi Multiome ATAC + Gene Expression dataset analyzed using Cell Ranger ARC 2.0.0

Mapped
95.7%

Assess data quality

View summary metrics to assess data quality and more.

View summary



Visualize and explore data

Discover differentially expressed genes, visualize your favorite genes, and explore your data with our visualization software.

Explore data

Dataset overview

Output and supplemental files


Input files

[Learn about Chromium analysis](#)

Batch download


Input files	File type	Size	md5sum
Sequencing data (FASTQ)	TAR	107 GB	0c5e9fb5964c38280361792cd444f2c6
Library (CSV)	CSV	407 B	b46af66895d55ca374a6376ce8fa1465

Step1. Generate single cell feature matrix




- **Alignment:** Mapping reads to the reference genome.
- **Peak Calling:** Identifying accessible chromatin regions
- **Fragment Counting:** Quantifying accessibility per region

Step1-1. Alignment




- **Alignment:** Mapping reads to the reference genome.
- **Peak Calling:** Identifying accessible chromatin regions
- **Fragment Counting:** Quantifying accessibility per region

Step1-2. Peak calling & Fragment counting



- Alignment: Mapping reads to the reference genome.
- **Peak Calling:** Identifying accessible chromatin regions
- **Fragment Counting:** Quantifying accessibility per region

Step1. Generate single cell feature matrix




Cell Ranger

Analyze Single Cell Gene Expression and Single Cell Immune Profiling data with a set of free, easy-to-use analysis pipelines.

[Learn more >](#) [Download >](#)

single cell RNA,
single cell Immune




Cell Ranger ATAC

Analyze Single Cell ATAC data with a set of free, easy-to-use analysis pipelines for identification of open chromatin regions, motif annotation, and differential accessibility analysis.

[Learn more >](#) [Download >](#)

single cell ATAC



Cell Ranger ARC

Analyze Single Cell Multiome ATAC + Gene Expression data with a set of free, easy-to-use pipelines for primary and secondary analyses.

[Learn more >](#) [Download >](#)

single cell Multiome
(ATAC +Gene expression)

<https://www.10xgenomics.com/software>

Step1-1. Download Cell Ranger ARC

<https://www.10xgenomics.com/support/software/cell-ranger-arc/downloads>

Runs on Linux systems

Cell Ranger ARC 2.0.2 (Aug 18, 2022)

- Chromium Single Cell Software Suite
- Self-contained, relocatable tar file. Does not require centralized installation.
- Contains binaries pre-compiled for CentOS/RedHat 7.0 and Ubuntu 14.04.
- Runs on Linux systems that meets the minimum compute requirements.

tar.gz compression

Download for Linux 64-bit (tar.gz)

File size: 699 MB

md5sum: 7303f8ceee7b60113c9a0087268830cd

curl wget

 Copy

```
curl -o cellranger-arc-2.0.2.tar.gz "https://cf.10xgenomics.com/releases/cell-arc/cellranger-arc-2.0.2.tar.gz"
```

References

[Steps to build references >](#)

[2020-A references & reference release notes >](#)

Human reference (GRCh38) - 2024-A

- Human reference (GRCh38) dataset required for Cell Ranger ARC.

Download Human Reference

File size: 14 GB


md5sum: 57a41ecf38b1ec2ef66b3d345ad05839

curl wget

 Copy

```
curl -O "https://cf.10xgenomics.com/supp/cell-arc/refdata-cellranger-arc-GRCh38-2024-A.tar.gz"
```


Step1-2. Specifying input FASTQ files for cellranger-arc count



Step1-2. Specifying input FASTQ files for cellranger-arc count

ATAC FASTQs

[Sample Name] *S1_L00* [Lane Number] [Read Type] _001.fastq.gz

Where **Read Type** is one of:

- **I1** : Dual index i7 read (optional)
- **R1** : Read 1
- **R2** : Dual index i5 read
- **R3** : Read 2

Cell Ranger ARC will also accept ATAC FASTQs in this format:

- **I1** : Dual index i7 read (optional)
- **R1** : Read 1
- **I2** : Dual index i5 read
- **R2** : Read 2

GEX FASTQs

[Sample Name] *S1_L00* [Lane Number] [Read Type] _001.fastq.gz

Where **Read Type** is one of:

- **I1** : Dual index i7 read (optional)
- **I2** : Dual index i5 read (optional)
- **R1** : Read 1
- **R2** : Read 2

Step1-3. Create a libraries CSV file

- CSV file format

Column Name	Description
fastqs	A fully qualified path to the <i>directory</i> containing the demultiplexed FASTQ files for this sample. This field does not accept comma-delimited paths. If you have multiple sets of fastqs for this library, add an additional row, and use the use same library_type value.
sample	Sample name assigned as the Sample_ID in the demultiplexing sample sheet.
library_type	This field is case-sensitive and must exactly match Chromatin Accessibility for a Multiome ATAC library and Gene Expression for a Multiome GEX library.

- example of libraries CSV file

```
fastqs,sample,library_type
/home/jdoe/runs/HNGEXSQXXX/outs/fastq_path,example,Gene Expression
/home/jdoe/runs/HNATACSQXX/outs/fastq_path,example,Chromatin Accessibility
```

<https://www.10xgenomics.com/support/software/cell-ranger-arc/latest/analysis/running-pipelines/single-library-analysis>

Step1-4. Run cellranger-arc count

cellranger-arc count --id=[A unique run ID string] \

--reference=[Path to the cellranger-arc-compatible reference package] \


--libraries=[Path to a 30column CSV file]

Output files


peak x count matrix

	Cell1	Cell2	...	CellN
Peak1	3	1	.	2
Peak2	2	0	.	1
Peak3	1	1	.	0
...
PeakM	5	1	.	0

fragment alignment file



QC info + QC summary in html



High-Quality Sample

Sample requiring attention due to suboptimal QC metrics

Due to the high likelihood of doublet formation, rigorous doublet filtering is recommended.

Sample A

Sample B



Alerts

The analysis detected 1 error.

Alert	Value	Detail
Number of cells detected is high	15,469	Ideally < 10,000 cells. This can be caused by incorrect quantification of the nuclei suspension, improper handling of nuclei,excessive background RNA and DNA, or unexpected behavior in the cell calling algorithm.




High-Quality Sample

Sample requiring attention due to suboptimal QC metrics


Cells ?

Estimated number of cells	6,913
Mean raw read pairs per cell	22,159.10
Fraction of high-quality fragments in cells	68.0%
Fraction of transposition events in peaks in cells	33.3%
Median high-quality fragments per cell	3,428



Cells ?

Estimated number of cells	8,531
Mean raw read pairs per cell	26,011.97
Fraction of high-quality fragments in cells	43.3%
Fraction of transposition events in peaks in cells	18.8%
Median high-quality fragments per cell	4,559




Alert	Value	Detail
⚠️ ATAC Transposition events within peaks is low	18.8%	It is expected that more than 25% of the transposition events fall within peak regions. A lower value could suggest a problem during library preparation causing excessive background transposition, low sequencing depth, or unexpected behavior in the peak calling algorithm resulting in very few peaks detected.
⚠️ ATAC Fragments in peaks is low	20.3%	Ideal > 25%. A low value can be caused by a problem during the transposition step, due to a population of cells with un-compacted DNA (e.g., activated granulocytes, dead or dying cells, unsupported organisms) or due to unexpected behavior in the peak calling algorithm resulting in very few peaks called.

High-Quality Sample

Targeting ?

Number of peaks	93,688
Fraction of genome in peaks	2.6%
TSS enrichment score	6.50
Fraction of high-quality fragments overlapping TSS	29.4%
Fraction of high-quality fragments overlapping peaks	36.0%



Sample requiring attention due to suboptimal QC metrics

⚠️ ATAC TSS enrichment 4.67 Ideal > 5. A low TSS score can be caused by poor sample prep, poor sample quality, a population of cells with highly accessible DNA (e.g., activated granulocytes, dead or dying cells, unsupported organisms) or poor quality reference genome annotation.

Targeting ?

Number of peaks	55,882
Fraction of genome in peaks	1.6%
TSS enrichment score	4.67
Fraction of high-quality fragments overlapping TSS	21.8%
Fraction of high-quality fragments overlapping peaks	20.3%

